Analytical Balance ES-J series

User Manual



Menu

One: Notes	1
Two: ES-J series balance features	1
Three: Balance appearance photo	2
Four: ES-J series balance technical parameters	
Five: ES-J series balance display and keypad drawings	4
Six: Installation	5
Seven: Function setting	6
Eight: Usage	10
Nine: Calibration	12
Ten: RS232 communication	12
Eleven: Attention of cleanliness	14
Twelve: Troubleshooting	15
Thirteen: Maintenance	15
Fourteen: Packing list	16

One: Notes

- 1.Please preheat the balance for 30 minutes in order to ensure the weighing is more accurate.
- 2. The normal work of the balance should be well-adapted to the environment, it should be put on the stable horizontal table. Don't put the balance in the unstable desk; do not put it at the place with a lot of dust in the air; do not put it in the vicinity of air-conditioning or blowing in the wind; do not put it in place of direct exposure under the sun; do not put it in place of heat flow; and do not put it in to generate and electric or magnetic field.
- 3.Please handle the heavy sample lightly while weighing, do not impact the weighing pan. If impacts it seriously, may lead to machenical system of the balance can not reset.
- 4.Please be careful while weighing liquid, avoid any liquid to inflow the internal balance from the bottom of weighing pan. If a similar situation occurs, should unplug the power supply immediately, clean internal liquid, or wait for the liquid to evaporate, ensure no residual liquid, then continue to use the balance.
- 5.A piece of lightly damp cloth is allowed to cover the balance after using.6.Please avoid any liquid on the surface of the balance.

Two: ES-J series balance features

- LCD liquid crystal display screen, more clear and bright, convenient for customers to read weighing value.
- All transparent glass windshield, samples can be seen entirely and clearly.
- Standard RS232 communication port—realize the communication with computer, printer or other device.
- Multiple weighing units change.

Three: Balance appearance photo



Four: ES-J series balance technical parameters

Test conditions: Room temperature 20°C, no air flow and no dust closed environment.

Model No.	ES-J120	ES-J220	ES-J120A	ES-J220A
Calibration	External		Internal	
Max capacity (g)	120	220	120	220
Readability (mg)		0	.1	
Repeatability (mg)		±().1	
Linearity error (mg)	±0.2			
Operational	12.25			
temperature range °C	13~25			
Operational moisture	10%~70%			
range RH	10%~/0%			
Response time	2.5			
(Average value) (s)	2.5			
Pan diameter (mm)	90			

External dimension (L*W*H)mm	230*310*330
Preheat time (minutes)	20-30
Baud rate options	300 600 1200 2400 4800 9600

Five: ES-J series balance display and keypad drawings

No.	Display information	Notes
1	g	Reading unit is "gram"
2	Ct	Reading unit is "carat"
3	Lb	Reading unit is "pound"
4	Oz	Reading unit is "ounce"
5	PCS	Count weighing state
6	%	Percentage weighing state
7	О	Data stable identification
8	CAL IN	Enter into calibration
9	CAL DN	Load calibration weight 100g
10	CAL	Calibrating
11	CAL-UP	Remove calibration weight
12	CAL-END	Calibration finish

13	Е	Overload max.capacity
14	-E	No weighing pan or incorrect weighing pan
15		"SET" key effective hint
16	•	Data processing
17	SAVE	Save
18	SEND	Save finish

Six: Installation

- 1.Installation of weighing pan's components: Install pan support and weighing pan accurately on the weighing pan support.
- 2.Adjust level: Adjust the two adjustable feet on the front bottom balance and be sure the level bubble in the middle part.
- 3. Plug the balance: Insert one end of power adapter into the 110-220V/50HZ AC plug base, and the other end of power adapter is inserted into the power plug base at rear part of the balance.
 - 4.Power on: After plugging the balance, model number "ESJ-4" will be displayed, 3 seconds later, display becomes black and displays power off state. Press ON/OFF key, "0.0000g" displays. If "P FAIL" displays, it means the balance was the previous non-normal shutdown, that is to shut down the balance without pressing OFF key. Press ON/OFF key, the

balance could be power on normally, displays "0.0000g".

Seven: Function setting

1 Back-light switch:

In the state of power-on, press SET key twice to open back-light. If the current state is open, that is to close back-light.

2 Parameter setting:

In the state of power-on, press SET key and then press PRINT key to display Cx-y, x is number 1-9, press TARE key to check circularly; y is the value set by Cx parameter, press PRINT key to change x. After setting, press OFF key, balance will display "SAVE……", press TARE key and "S-END" displays, wait for about 2 seconds, save finishes and updates data and back to weighing state. Press ON/OFF key twice to take effect new parameter. Please refer to below list as details:

Cx	Сх-у	Description
	C1-0	Auto-calibration (just for internal calibration balance)
C1. C-13	C1-1	Manual external weight calibration
C1: Calibration way	C1-2	Manual internal calibration
	C1-3	Adjustment between internal and external calibration
C2: Basic sample	C2-0	Quantity is 10

data option in count	C2-1	Quantity is 20
weighting state	C2-2	Quantity is 50
	C2-3	Quantity is 100
	C2-4	Quantity is 1000
	C3-0	No zero tracking
	C3-1	1d
	C3-2	2d
C3: Zero tracking	C3-3	3d
	C3-4	4d
	C3-5	5d
	C3-6	Can not be operated by user
	C4-0	1200bps
C4: Serial port baud	C4-1	2400bps
rate option	C4-2	4800bps
	C4-3	9600bps
	C5-0	Return to zero stable output
C5: Data output way	C5-1	Stable output
option	C5-2	Command control output
	C5-3	Continuous output

	C5-4	Press PRINT key output
C4. Vov. sound	C6-0	no
C6: Key sound	C6-1	yes
	C7-0	Low
C7:Anti-interference	C7-1	Medium
degree	C7-2	High
	C7-3	Strong

* Return to zero stable output:

Load weight when displays "0.0000g", after stable, output stable data to external device.

* Stable output:

When data is stable again, output current stable data.

* Command control output:

The P<CR<LF> commands which are outputed by the PRINT key on the balance through serial port are regarded as "print control command". There will be a group of current weight data output after the balance receives a command.

* Continuous output:

The balance is fixed to send the weight data every 0.3 second.

3 Command control:

Effective commands are as follows:

①O<CR><LF> ON/OFF command: same function with ON/OFF key on the balance (4F 0D 0A);

②T<CR><LF> TARE command: same function with TARE key on the balance (54D 0D 0A);

③C<CR><LF> mode conversion command: same function with CAL key on the balance (43 0D 0A);

④M<CR><LF> mode conversion command: same function with MODE key on the balance (4D 0D 0A);

⑤P<CR><LF> print command: same function with PRINT key on the balance. A group of current weight data will be output to external device as soon as the balance receives this command (50D 0D 0A);

<CR><LF> descriptions are as follows:

<CR>: Return symbol (0D)

<LF>: Enter symbol (0A)

Eight: Usage

1. Ordinary weighing:

(1) After plugging, in order to achive best weighing effect, balance should be preheated at least 30 minutes.

- (2) Keep weighing pan clean, press ON/OFF key and "0.0000g" displays.
 - (3) If other weighing units or weighing modes are needed, press MODE key to adjust the displayed data to other units or weighing mode data.
 - (4) Open the door of weighing chamber, put sample onto the middle of weighing pan gently, then close the door and wait for reading data after stable.
- (5) Open the door of weighing chamber to remove sample for the next weighing; if do not weigh, please close the door to avoid dust into the balance.

2. Container weighing:

- (1) Put the container on the weighing pan;
- (2) Press TARE key after "O" (stable mark) displays, "0.0000g" is displayed;
 - (3) Put sample into the container;
 - (4) Read the value of sample after "O" (stable mark) displays;

3. Count weighing:

- (1) Choose sample's quantity according to system parameter list;
- (2) Press TARE key and "0.0000g" displays after stable;
- (3) Press MODE key to "Count Mode" state;

- (4) Put sample in the middle of weighing pand and close the door of weighing chamber;
 - (5) Press SET key and balance system will take sample according to C2 parameter;
- (6) After sampling, balance will display the sample's quantity according to C2 parameter. Remove the sample and wait for balance to stable zero, user can operate count weighing.

Attention: The minimum readable value of sample quantity could not be less than the minimum resolution of balance.

4. Percentage weighing:

- (1) Press TARE key and "0.0000g" displays after stable;
- (2) Press MODE key to "Percentage Mode" state;
- (3) Put sample in the middle of weighing pan and close the door of weighing chamber;
- (4) Press SET key and balance system will take this reference as 100.00% basic value;
- (5) After sampling successfully, "100%" displays. Remove sample and after balance returns back to stable zero, user can operate percentage weighing.

Attention: The minimum readable value of sample quantity could not be

less than the minimum resolution of balance.

Nine: Calibration

1. External calibration

Menu setting (C1-1)

100g External weight calibration:

(1) Remove sample out of weighing pan, close the door of weighing

chamber, "0.0000g" displays;

(2) Press CAL key and "CAL IN" – "CAL …" displays, finally "CAL

dn" displays several seconds later;

(3) Put F1 grade 100g standard weight on the middle of weighing pan

and close the door of weighing chamber;

(4) When "CAL UP" displays (5-8 seconds later), remove weight;

(5) Close the door of weighing chamber and wait for 5-8 seconds,

"CAL-END" displays and calibration is finished;

(6) If the deviation is more than ± 1 between displayed weight value

and calibration weight, re-calibration is needed till the deviation is within

 ± 1 , that means calibration is finished successfully.

200g External weight calibration:

(1) Remove sample out of weighing pan, close the door of weighing

11

chamber, "0.0000g" displays;

- (2) Put F1 grade 100g standard weight on the middle of weighing pan and close the door of weighing chamber;
- (3) Press SET key and "-----" displays, then press MODE key and "CAL 200" displays;
- (4) The calibrated value will be displayed after 5-8 seconds, calibration is finished.

2. Automatic internal calibration (just for internal calibration balance)

Menu setting (C1-0).

Automatic adjustment is performed:

- (1) After the balance is turned on, it will undergo a fully automatic calibration.
- (2) Having detected a change in temperature by 2°C from the previous adjustment,
- (3) A defined time interval has elapsed, selectable 2, 3, 4, 6, 8 hours (C7-0 C7-4).

Attention: During the calibration process of the balance, it is prohibited to perform any operations on the balance.

In case of the adjustment error (e.g. object on the plate) the following

error message will be displayed: <CAL NO>, manually repeat the adjustment procedure.

- 3. Manual internal calibration (just for internal calibration balance) Menu setting (C1-0 / C1-2)
- (1) Press CAL button in the weighing mode.
- (2) The adjustment is carried out automatically.
- (3) After the successful adjustment the scales will be automatically switched back to the weighing mode.

In case of the adjustment error (e.g. object on the plate) the following error message will be displayed: <CAL NO>, repeat the adjustment procedure.

4. Adjustment between internal and external calibration

Menu setting (C1-3)

- (1) In the weighing mode, press the CAL button, the balance displays "CAL 100", then the balance shows "", followed by "CAL dn".
- (2) Then place the standard E2 weight of 100g on the balance pan, the balance displays "" indicating that it is comparing, and then the balance shows "CAL up". At this point, remove the E2 weight of 100g, the balance displays "______", and then shows "-----", the calibration is complete and it returns to the weighing mode.

Note: Before calibrating the balance, please adjust the level first to ensure

the balance is in a horizontal state. Allow the balance to warm up for more

than 30 minutes after powering on.

Ten: RS232 communication

1. Serial port working mode: Induction communication

Baud rate: please refer to system parameter list C4 parameter

Start bit: 1

Data bit: 8

Check bit: no

Stop bit: 1

2. Weight data transmitting mode: please refer to system parameter list C5

parameter

3. Weight data transmitting format: FXXXX.XXXXKKK<CR><LF> have

the following significance:

F: sign bit "+" or "-"

X: 0-9 weight data

.: decimal point

K: reserved three bit weighing unit symbol, right side align, spaces

added if less than three bits

<CR>: Return symbol

14

<LF>: Enter symbol

If weight data is +10.0000g, the data format of serial port is: +0010.0000g<CR>LF>

4.Decision method of stable and unstable data output:

When unstable data outputs, non-unit indicating bit of data serial port,

KKK three bits are spaces; when stable data outputs, unit information will
be output.

5.External serial data line connection method:

Electronic	Computer	Electronic	Serial port
Balance	Computer	Balance	Printer
9 needles	 9 holes	9 needles	 25 needles
2	 2	2	 2
3	 3	3	 3
5	 5	5	 7

Connection diagram of electronic balance and computer, printer

Attention: length of data line should not exceed 15 meters.

Eleven: Attention of cleanliness

- Unplug power line before cleaning;
- Suggest use alcohol or soft solvents instead of corrosive cleaning

detergent;

- •Do not let water or other liquid splashed into the balance while cleaning;
- Use a dry soft cloth to dry off the balance after cleaning.

Twelve: Troubleshooting

Failure phenomenons	Failure reasons	Solutions
Can not power on	•Do not connect power supply •AC/DC damaged	•Plug power cord •Change power adapter
	•Working environment is bad	•Improve working environment to avoid vibration and air disturbance
	•Door of weighing chamber is not closed tightly	• Close the door tightly
Displayed value is unstable • Foreign object or scratch between weighing pan and housing of balance • Power supply is Unstable and exceed allowable value		• Remove foreign object, turn weighing pan to prevent scratch • Access 220V AC
	• Measured sample is unstable (e.g. absorb water or evaporate)	power regulator
Displayed value is	Balance is not calibrated	Calibration
out of accord with	Do not tare before	• Press TARE key

actual weight	weighing	
value	Do not adjust level	Adjust level feet

Thirteen: Maintenance

We will guarantee the balance within one year from the date of sale except the following status:

- 1. The guarantee period is expired;
- 2. The balance is damaged due to user's own negligence.
- 3. The balance is damaged due to the exposure to the environment which has radioactive or corrosive material.

Fourteen: Packing list

No.	Description	Quantity
1	Balance main body	1
2	Weighing pan	1
3	Weighing pan support	1
4	English user manual	1
-	Standard weight (F1	
5	grade) 100g	1